



Anti-SH3GL1 SH3-domain Grb2-like protein, NeuroMab clone K51/1

Immunogen:

Synthetic peptide amino acids 287-302 of rat SH3GL1 (accession number AAH70893)

Human: 81% identity (13/16 amino acids)

Mouse=81% identity (13/16 amino acids)

Monoclonal antibody info:

Mouse strain: Balb/C

Myeloma cell: SP2/0

Mouse Ig Isotype: IgG1

NeuroMab Applications:

Immunoblot, Immunohistochemistry and Immunoprecipitation

Species Reactivity: human, mouse, rat

MW: 43 kDa

Adult rat brain membrane immunoblot

